

14ml


Round Bottom Tube

17 x 95 mm style

This is made from Polystyrene high grade material.

- High Strength,
- High Resistance to Chemical,
- High Resistance to Impact
- High Resistance to Heat.

It has very clear and transparent surface, and MEA certified ensuring Compliance in IVF application. It has been produced in fully automated high quality finish mould. It has been produced in absolutely clean room and it is certified RNase, DNase, Human DNA and pyrogen free.



**MEA
CERTIFIED**
one cell MEA \geq 80%
Expanded Blastocyst
within 96H
Embryotech
- USA



Cat No	QTY / Case	Description
400020004	250 nos	14ml Round Bottom Tubes-17X95mm Style With Snap Cap-Individually Packed Sterile By Gamma radiation & MEA Certified
400020008	300 nos	14ml Round Bottom Tubes-17X95mm Style With Snap Cap-(Pack of 5 Sterile) By Gamma radiation & MEA Certified

It can be easily withstand centrifugal force up to 4000 x g.



APPLICATIONS

Cell separation | IUI | IVF | Sample storage | Liquid handling



**ARCK
TECHNOLOGIES**
INNOVATIONS PVT LTD

📍 **Manufactured By :**
Arck Technologies Innovation Pvt. Ltd.
B2 - 704, Boomerang, Chandivali farm Road,
Powai, Mumbai - 400072

✉ **E mail :** info@arcktechnologies.com

☎ **Tel no :** +91 77770 89451, +91 84339 68529