

# 15ml Graduated Conical Centrifuge Tube (17mm X 120 mm)



- 1) It is made of high clarity, high grade Polypropylene.
- 2) It is Manufactured in RNase, DNase and Pyrogen free room.
- 3) It has very clear transparent Surface.
- 4) It can be centrifuged up to 16000 RCF.
- 5) It has leak proof screw cap.
- 6) It has large writing area of sample identification.
- 7) It is Sterile.



Cat No	QTY / Case	Description
40002005	250 nos	15ml Graduated Centrifuge Tubes-17X120mm Sterile-Polypropylene

## APPLICATIONS

Cell separation | IUI | IVF | Sample storage | Liquid handling



**ARCK**  
**TECHNOLOGIES**  
INNOVATIONS PVT LTD

📍 **Manufactured By :**  
Arck Technologies Innovation Pvt. Ltd.  
B2 - 704, Boomerang, Chandivali farm Road,  
Powai, Mumbai - 400072

✉ **E mail :** [info@arcktechnologies.com](mailto:info@arcktechnologies.com)

☎ **Tel no :** +91 77770 89451, +91 84339 68529