

# 3ML PASTEUR (TRANSFER) PIPETTES



**ARCK  
TECHNOLOGIES**  
INNOVATIONS PVT LTD



3 ml Pasteur pipettes are made from Low Density Polyethylene (LDPE) Hydrophobic Material.

The pipettes are clear, translucent and 6" inch (15.4 cm ) in length.

Produced in fully automated high-quality finish mould and clean room production ensures product is free from RNase, DNase, Endotoxin free.




Gamma radiation sterilization ensures optimal sterility and MEA certified ensures compliance for IVF application.

Graduated at 0.5 ml, 1.00 ml, 1.5 ml, 2.00 ml, 2.5 ml and 3.00 ml and can be easily visible.

The tip has been designed in such a way that samples are safe during the procedure, and shows consistent reproduction of drop size.

These high-quality pipettes ideal for a wide variety of sample preparation, media mixing etc.

## APPLICATIONS

-  IVF
-  IUI
-  Liquid handling



Cat No	Description	Pack Size
TR003001	3 ml Pasteur (Transfer) Pipette <b>Individual packed, Gamma Sterile</b> 6' inch (15.4 cm) <b>MEA Certified from -Embryotech -USA</b>	300 nos



**Manufactured By :**  
Arck Technologies Innovation Pvt. Ltd.  
B2 - 704, Boomerang, Chandivali farm Road,  
Powai, Mumbai - 400072

✉ E mail : [info@arcktechnologies.com](mailto:info@arcktechnologies.com)  
☎ Tel no : +91 77770 89451, +91 84339 68529