



ARCK 4-Well Dish


The 4-Well Dish is mainly used for oocyte & embryo handling & culture

It is made from Medical Grade Polystyrene material & certified for FDA-21 CFR

Our 4-well dish is designed & developed with innovative features for quality & ease of use.

- Ensures proper well orientation & prevents cross contamination with notched lid & dish.
- It is MEA Certified & Endotoxin-LAL Certified from Well Known certified lab in USA.
- True flat bottom for maintaining optimal temperature.
- Easy grip design for increased safety.
- Available in Individual sterile pack
- It is manufactured in complete clean room ensuring it is free from Human DNA, RNA Free & Endotoxin Free.
- It is Gamma radiation sterilized to ensure sterility of the product.

Cat No	QTY / Case	Description
40002007	250 nos	ARCK 4- well Dish, 1.9cm ² /Well, Individually Pack, Sterile (MEA Certified)



ARCK
TECHNOLOGIES
INNOVATIONS PVT LTD

Manufactured By : Arck Technologies Innovation Pvt. Ltd.
B2 - 704, Boomerang, Chandivali farm Road, Powai, Mumbai - 400072

✉ E mail : info@arcktechnologies.com
☎ Tel no : +91 77770 89451, +91 84339 68529