

5ml Round Bottom Tube

12 x 75 mm style



**ME A
CERTIFIED**
one cell MEA \geq 80%
Expanded Blastocyst
within 96H
Embryotech
- USA



This is made from Polystyrene high grade material.

- High Strength,
- High Resistance to Chemical,
- High Resistance to Impact
- High Resistance to Heat.

It has very clear and transparent surface, and MEA certified ensuring Compliance in IVF application. It has been produced in fully automated high quality finish mould. It has been produced in absolutely clean room and it is certified RNase, DNase, Human DNA and pyrogen free.

Cat No	40002003	Description
Pack Size	Individually packed	5ml Round Bottom Tubes-12X75mm Style With Snap Cap- Individually Packed Sterile By Gamma radiation & MEA Certified
QTY / Case	500 nos	

It can be easily withstand centrifugal force up to 4000 x g.



APPLICATIONS

Cell separation | IUI | IVF | Sample storage | Liquid handling

📍 **Manufactured By :**
Arck Technologies Innovation Pvt. Ltd.
B2 - 704, Boomerang, Chandivali farm Road,
Powai, Mumbai - 400072

✉ **E mail :** info@arcktechnologies.com

☎ **Tel no :** +91 77770 89451, +91 84339 68529



**ARCK
TECHNOLOGIES**
INNOVATIONS PVT LTD