

# Culture Dish-Sterile (90mm x 20mm)


- 90mm culture dish is mainly used for scanning for follicular fluid
- It can be used for Embryo & Gamete culture
- Oocyte insemination & IVF Denudation



**It is made from Medical Grade Polystyrene material & certified for FDA-21 CFR**

Our 90mm dish is designed & developed with innovative features for quality & ease of use.

- 1) True flat bottom for maintaining optimal temperature.
- 2) Easy grip design for increased safety.
- 3) Available in pack of 5.
- 4) Culture dish with Airvent making sure proper gas exchange takes place & also prevent contamination.
- 5) It is manufactured in complete clean room ensuring it is free from Human DNA, RNA Free & Endotoxin Free.
- 6) It is Gamma radiation sterilized to ensure sterility of the product.
- 7) It is MEA Certified & Endotoxin-LAL Certified from Well Known certified lab in USA.



Cat No	QTY / Case	Description
40002016	330 nos	Culture/Petri Dish 90mm x 20mm, Pack of 5, Sterile (MEA Certified)



**ARCK  
TECHNOLOGIES**  
INNOVATIONS PVT LTD

📍 **Manufactured By :**  
Arck Technologies Innovation Pvt. Ltd.  
B2 - 704, Boomerang, Chandivali farm Road,  
Powai, Mumbai - 400072

✉ **E mail :** [info@arcktechnologies.com](mailto:info@arcktechnologies.com)

☎ **Tel no :** +91 77770 89451, +91 84339 68529